

# NICORETTE® Inhaler

Reduces cravings and satisfies the hand-to-mouth ritual

## How to properly use the NICORETTE® Inhaler

1

Remove the mouthpiece from the plastic wrap. Align marks to separate mouthpiece in two.




2

Peel back foil and remove one cartridge from tray. Press firmly into bottom of mouthpiece until the seal breaks.



3

Replace top of mouthpiece. Align marks and press down firmly to break the top foil seal of cartridge. Twist and secure.



4

Place the tapered end of the inhaler in your mouth and inhale deeply into the back of your throat or puff in short breaths.

*To ensure safe disposal, used cartridges should be discarded out of reach children and pets.*

Patients should not concurrently smoke or use any other nicotine-containing products. The Inhaler should not be used for more than 20 consecutive minutes per session.

© Johnson & Johnson Inc. 2015.